Next-FACS planing-upon the insights on the GAL1 promoter

Contents

LIUIC .	14012020-	INCA	U-I	AC	ս I	Jia	TIII	11E	,- u	ιP	911	. 01	iic	11	191	81	103	, 0	11	UI.	ıc	u	 	ŀ	,,,	<i>J</i> 11	ııc	,,,	, I			
Date	e																														 	
Obje	ective																										. .				 	
Met	hod																															
Resu	ılts																															
Con	clusion																															
List	of Figui	res																														
1	Experiment	tal De	esig	n fo:	r tl	he	mi	icro	osc	or	ΟV	me	eas	sur	en	ıeı	$_{ m nts}$, ,														

1

Title: 14012020 Novt-FACS planing-upon the insights on the CAL1 promotor

Title: 14012020- Next-FACS planing-upon the insights on the GAL1 promoter

Date

14012020

Next experiment date:

Objective

To set up a consistent FACs experiment with the microscopy experiments I did with Werner strains containing the Gal1p-Cdc42-sfGFP.

Method

- Insights from the Gal1 promoter
 - The location in the genome at which the gal promoter was inserted can have a strong effect on the expression pattern of the gal promoter (*Ramon feedback*). Hence we should not compare different studies of the Gal1 promoter with ours if the integration in the genome is in a different location and also if it is a plasmid or not.
 - We should compare systematically the Gal1p expression pattern of the strains that has the sfGFP (Werner strains) and the mneonGreen ones (Ramon/Miranda strains), because they have the same type of genomic integration of the Gal1 promoter.
 - * Ask Reza for his data with WT+mneongreen to comapre with mine
- Follow the same protocol as I followed for the microscopy measurements




Figure 1: Experimental Design for the microscopy measurements

- 15 hours of incubation in 2% Gal +2% Raff (i.e. overnight incubation from 17:00 to 08:00)
- Washing step with CSM+2% Raff+0%Gal to the respective Galactose (at 08:00-09:00) concentrations. Incubate
- Measure FACs after 24 hours of incubation. (next day of the washing step)
 - Use the references cdc42-GFP ywkd038 and ywkd001 in 2% dextrose +2% Raff, in CSM-met and CSM respectively.

Results

Conclusion